

PCR protocol

RBRC 04360

Strain name C57BL/6J-Tg(hSMA-TRPC6)1

Primer	primer_name	primer_seq	length	
1	FHSA-Bam	TGCTCCAGCCTTTCTAATTTTATG	24	Depositor's primer
2	RTrpc6-PCR	CTGGTCCTCGATTAGCTAACCTTCT	25	Depositor's primer
3	oIMR0042	CTAGGCCACAGAATTGAAAGATCT	24	JAX Genotyping Protocol
4	oIMR0043	GTAGGTGGAAATTCTAGCATCATCC	25	JAX Genotyping Protocol


Reaction Components

Enzyme TaKaRa Ex Taq (Cat. No. RR001A)

Reaction Conditions

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		11.3			
PCR buffer	10 X	2.0			
dNTP	2.5 mM	1.6			
primer 1	10 pmol/ μ l	1.0			
primer 2	10 pmol/ μ l	1.0			
primer 3	10 pmol/ μ l	1.0			
primer 4	10 pmol/ μ l	1.0			
DNA		1.0			
Taq		0.1			
total		20			

step	temp (°C)	time (sec)	note
1	94	180	
2	94	30*	
3	60	30*	
4	72	90*	*Repeat for 30 cycles
5	72	300	



Marker

1 kb DNA ladder (Invitrogen, Cat. 15615-024)

Gel

1% agarose in 1 X TAE (100 V, 30 min)

PCR products

primer_set		product_size	product
primer 1	primer 2	1300 bp	Tg
primer 3	primer 4	324 bp	Internal control (IL-2)

Comment